

PRODUCT DATA SHEET

Anti-RAB5 Rabbit Polyclonal antibody

Cat Code: HBM-RAB5-PR1-###

Size: 100 or 50 µg

Description	Rab5 is a small GTPase belonging to Ras superfamily of monomeric G proteins. Rab GTPases play an essential role in the regulation of membrane traffic and are involved in vesicle formation and trasport and fusion to the membrane. Rab5 is expressed on exosome membranes and it might have an active role during endo/exocytotic processes of microvesiscles through the plasma membrane
Immunogen	Full length of RAB5 protein
Туре	Rabbit Polyclonal
Species Reactivity	Human.
Concentration	0,2 mg/ml.
Application	ELISA, Flow Cytometry, Western blotting
Recommended dilution	Recommended dilution for ELISA and WB 1/500
Storage buffer	Glycine 0.1 M, NaCl 0.5 M, Tris-HCl 0.1 M with 15 mM sodium azide, pH 7.4
Storage Condition	Store at 4°C for short term . Over 1 month store at -20°C . Avoid freezing and thawing cycles.
Expiration	See vial label
Lot Number	See vial label
Product notes	Sodium azide yields highly toxic hydrazoic acid under acidic conditions. Dilute azide compounds in running water before discarding to avoid accumulation of potentially explosive deposits in plumbing.

This product is for research use only and is not approved for use in humans or in clinical diagnosis.

HansaBioMed LLC Akadeemia tee 15, 12618 Tallinn, ESTONIA

www.hansabiomed.eu Tel: +372 6204348 Fax: +372 6204358 email: info@hansabiomed.eu



Contact in USA Galen Laboratory Supplies P.O. Box 682 Middletown, CT 06457. USA

www.galenlabsupplies.com Tel: 860-704-9058 Fax: 888-260-6091 email: info@ galenlabsupplies.com