



Anti-DNMT3A Antibody (Clone ABM1F35)

Alternative Names: DNMT3A, DNA (cytosine-5)-methyltransferase 3A, DNA methyltransferase HsallIA, DNA MTase HsallIA, M.HsallIA, Dnmt3a

Catalogue Number: AA17-10027-50ug

Size: 50 µg

Background Information

DNA methyltransferase 3A (DNMT3A) belongs to the mammalian methyltransferase gene family which is responsible for tissue-specific gene expression. DNMT3A together with other methyltransferases carries out de novo methylation of cytosine residues in CpG islands by the enzymatic addition of methyl residues from S-adenosyl-L-methionine to the 5-carbon atom of the cytosine ring.

DNMT3A also represses key genes enabling the coupling of insulin secretion to glucose levels during beta cell maturation. Loss of Dnmt3a in pancreatic beta cells prevents this developmental metabolic reprogramming, resulting in loss of GSIS (Glucose-Stimulated Insulin Secretion). Dnmt3a overexpression has been shown in several cancers including AML (Acute Myeloid Leukemia).

Product Information

Antibody Type:	Monoclonal	Host:	Mouse
Isotype:	IgG1 kappa	Species Reactivity:	Human
Immunogen:	Full length recombinant human DNMT3a		
Format:	50 µg in 100 µl PBS containing 0.05% BSA and 0.05% sodium azide.		
Storage Conditions:	6 months: 4°C. Long-term storage: -20°C. Avoid multiple freeze and thaw cycles.		
Applications:	WB IHC WB: 2-4 µg/ml, Immunohistochemical analysis: 5 µg/ml		

Additional Information

Subcellular location:	Nucleus	MW:	102kDa (Intended as a general guide and does not allow for all isoforms and species variations)
Gene ID	1788	Uniprot ID:	Q9Y6K1